

产品名称: GSK3 α/β Rabbit Monoclonal Antibody
产品货号: AMRe21392



产品概述 (Summary)

产品名称 (Production Name)	GSK3 α/β Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG,Kappa
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Protein A

免疫原信息 (Immunogen)

基因名 (Gene Name)	GSK3A/GSK3B
别名 (Alternative Names)	GSK3A;Glycogen synthase kinase-3 α ;GSK-3 α ;Serine/threonine-protein kinase GSK3A;GSK3B;Glycogen synthase kinase-3 β ;GSK-3 β ;Serine/threonine-protein kinase GSK3B
基因 ID (Gene ID)	2931
蛋白 ID (SwissProt ID)	P49840;P49841.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:2000-1:10000,IHC 1:400-1:800,ICC/IF 1:200-1:1000,ELISA 1:5000-1:20000,IP 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:51kD,47kD;Observed MW:51kD,47kD

产品名称: GSK3 α/β Rabbit Monoclonal Antibody
产品货号: AMRe21392




研究背景 (Background)

Cell localization: Cytoplasm. glycogen synthase kinase 3 alpha (GSK3A) Homo sapiens This gene encodes a multifunctional Ser/Thr protein kinase that is implicated in the control of several regulatory proteins including glycogen synthase, and transcription factors, such as JUN. It also plays a role in the WNT and PI3K signaling pathways, as well as regulates the production of beta-amyloid peptides associated with Alzheimer's disease. [provided by RefSeq, Oct 2011],

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from A549

cells, using GSK3 α/β Rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG antibody was used to detect the antibody.

注意事项 (Note)

For research use only .