

产品名称: Catenin δ -1 Rabbit Monoclonal Antibody
产品货号: AMRe21399



产品概述 (Summary)

产品名称 (Production Name)	Catenin δ -1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG,Kappa
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Protein A

免疫原信息 (Immunogen)

基因名 (Gene Name)	CTNND1 KIAA0384 Cadherin associated Src substrate;Cadherin-associated Src substrate;CAS;Catenin;cadherin associated protein;delta 1;Catenin delta 1;Catenin delta;Catenin delta-1;CTND1_HUMAN;CTNND 1;CTNND;CTNND1;delta 1 Catenin;KIAA0384;p120;P120 CAS;p120 catenin;P120 CTN;p120;cas;p120;ctn;P120CAS;P120CTN
别名 (Alternative Names)	
基因 ID (Gene ID)	1500.0
蛋白 ID (SwissProt ID)	O60716.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:200-1:1000,ICC/IF 1:200-1:1000,ELISA 1:5000-1:20000,IP 1:50-1:200
----------------------	--

产品名称: Catenin δ -1 Rabbit Monoclonal Antibody
产品货号: AMRe21399




蛋白分子量 (Molecular Weight) Calculated MW:108kD;Observed MW:108kD

研究背景 (Background)

Cell localization:Membrane.catenin delta 1(CTNND1) Homo sapiens This gene encodes a member of the Armadillo protein family, which function in adhesion between cells and signal transduction. Multiple translation initiation codons and alternative splicing result in many different isoforms being translated. Not all of the full-length nature of the described transcript variants have been determined. Read-through transcription also exists between this gene and the neighboring upstream thioredoxin-related transmembrane protein 2 (TMX2) gene. [provided by RefSeq, Dec 2010],

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from 3T3-L1

cells, using Catenin δ -1 Rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG antibody was used to detect the antibody.

注意事项 (Note)

For research use only .