

产品名称: MelanA Rabbit Monoclonal Antibody
产品货号: AMRe21411



产品概述 (Summary)

产品名称 (Production Name)	MelanA Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA,IP
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG,Kappa
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Protein A

免疫原信息 (Immunogen)

基因名 (Gene Name)	MLANA MART1
别名 (Alternative Names)	Melanoma antigen recognized by T-cells 1;MART-1;Antigen LB39-AA;Antigen SK29-AA;Protein Melan-A;
基因 ID (Gene ID)	2315.0
蛋白 ID (SwissProt ID)	Q16655.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:400-1:1000,ICC/IF 1:200-1:1000,ELISA 1:5000-1:20000,IP 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:13kD;Observed MW:18kD

研究背景 (Background)


产品名称: MelanA Rabbit Monoclonal Antibody
产品货号: AMRe21411



Cell localization:Endoplasmic reticulum membrane.tissue specificity:Expression is restricted to melanoma and melanocyte cell lines and retina.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from B16F10

cells, using MelanA Rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG antibody was used to detect the antibody.

注意事项 (Note)

For research use only .