

产品名称: MEK1/2 Rabbit Monoclonal Antibody
产品货号: AMRe21418



产品概述 (Summary)

产品名称 (Production Name)	MEK1/2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG,Kappa
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Protein A

免疫原信息 (Immunogen)

基因名 (Gene Name)	MAP2K1/MAP2K2
别名 (Alternative Names)	MAP2K1;MEK1;PRKMK1;Dual specificity mitogen-activated protein kinase kinase 1;MAP kinase kinase 1;MAPKK 1;MKK1;ERK activator kinase 1;MAPK/ERK kinase 1;MEK 1;MAP2K2;MEK2;MKK2;PRKMK2;Dual specificity mitogen-activated protein k
基因 ID (Gene ID)	5604;5605
蛋白 ID (SwissProt ID)	Q02750;P36507.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:500-1:2000,ICC/IF 1:200-1:1000,ELISA 1:5000-1:20000,IP 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:44kD;Observed MW:44kD

产品名称: MEK1/2 Rabbit Monoclonal Antibody
产品货号: AMRe21418




研究背景 (Background)

Cell localization: Cytoplasm. The protein encoded by this gene is a member of the dual specificity protein kinase family, which acts as a mitogen-activated protein (MAP) kinase kinase. MAP kinases, also known as extracellular signal-regulated kinases (ERKs), act as an integration point for multiple biochemical signals. This protein kinase lies upstream of MAP kinases and stimulates the enzymatic activity of MAP kinases upon wide variety of extra- and intracellular signals. As an essential component of MAP kinase signal transduction pathway, this kinase is involved in many cellular processes such as proliferation, differentiation, transcription regulation and development. [provided by RefSeq, Jul 2008],

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from Hela

cells, using MEK1/2 Rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG antibody was used to detect the antibody.

注意事项 (Note)

For research use only .