

产品名称: IKK α/β Rabbit Monoclonal Antibody
产品货号: AMRe21420



产品概述 (Summary)

产品名称 (Production Name)	IKK α/β Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG,Kappa
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Protein A

免疫原信息 (Immunogen)

基因名 (Gene Name)	CHUK/IKBKB CHUK;IKKA;TCF16;Inhibitor of nuclear factor kappa-B kinase subunit alpha;I-kappa-B kinase alpha;IKK-A;IKK-alpha;IkBKA;IkappaB kinase;Conserved helix-loop-helix ubiquitous kinase;I-kappa-B kinase 1;IKK1;Nuclear factor NF-kappa-B
别名 (Alternative Names)	
基因 ID (Gene ID)	1147;3551
蛋白 ID (SwissProt ID)	O15111;O14920.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:500-1:2000,ICC/IF 1:200-1:1000,ELISA 1:5000-1:20000,IP 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:86kD;Observed MW:86kD

产品名称: IKK α/β Rabbit Monoclonal Antibody
产品货号: AMRe21420




研究背景 (Background)

Cell localization: Cytoplasm, Nucleus. This gene encodes a member of the serine/threonine protein kinase family. The encoded protein, a component of a cytokine-activated protein complex that is an inhibitor of the essential transcription factor NF-kappa-B complex, phosphorylates sites that trigger the degradation of the inhibitor via the ubiquitination pathway, thereby activating the transcription factor. [provided by RefSeq, Jul 2008],

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from HEK293

cells, using IKK α/β Rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG antibody was used to detect the antibody.

注意事项 (Note)

For research use only .