

产品名称: Paxillin Rabbit Monoclonal Antibody
产品货号: AMRe21430



产品概述 (Summary)

产品名称 (Production Name)	Paxillin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG,Kappa
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Protein A

免疫原信息 (Immunogen)

基因名 (Gene Name)	PXN
别名 (Alternative Names)	PXN;Paxillin
基因 ID (Gene ID)	5829.0
蛋白 ID (SwissProt ID)	P49023.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:1000-1:4000,ICC/IF 1:200-1:1000,ELISA 1:5000-1:20000,IP 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:68kD;Observed MW:68kD

研究背景 (Background)

Cell localization:Cytoplasm.This gene encodes a cytoskeletal protein involved in actin-membrane attachment at sites of cell


产品名称: Paxillin Rabbit Monoclonal Antibody
产品货号: AMRe21430



adhesion to the extracellular matrix (focal adhesion). Alternatively spliced transcript variants encoding different isoforms have been described for this gene. These isoforms exhibit different expression pattern, and have different biochemical, as well as physiological properties (PMID:9054445). [provided by RefSeq, Aug 2011],

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from Hela

cells, using Paxillin Rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG antibody was used to detect the antibody.

注意事项 (Note)

For research use only .