

产品名称: Adiponectin Rabbit Monoclonal Antibody
产品货号: AMRe21479



产品概述 (Summary)

产品名称 (Production Name)	Adiponectin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG,Kappa
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Protein A

免疫原信息 (Immunogen)

基因名 (Gene Name)	ADIPOQ
别名 (Alternative Names)	ADIPOQ;ACDC;ACRP30;APM1;GBP28;Adiponectin;30 kDa adipocyte complement-related protein;Adipocyte complement-related 30 kDa protein;ACRP30;Adipocyte;C1q and collagen domain-containing protein;Adipose most abundant gene transcript 1
基因 ID (Gene ID)	9370.0
蛋白 ID (SwissProt ID)	Q15848.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:1000-1:4000,ICC/IF 1:200-1:1000,ELISA 1:5000-1:20000,IP 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:26kD;Observed MW:30kD

产品名称: Adiponectin Rabbit Monoclonal Antibody
产品货号: AMRe21479




研究背景 (Background)

Cell localization: Secreted. adiponectin, C1Q and collagen domain containing (ADIPOQ) Homo sapiens This gene is expressed in adipose tissue exclusively. It encodes a protein with similarity to collagens X and VIII and complement factor C1q. The encoded protein circulates in the plasma and is involved with metabolic and hormonal processes. Mutations in this gene are associated with adiponectin deficiency. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Apr 2010],

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from Rat fat

cells, using Adiponectin Rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG antibody was used to detect the antibody.

注意事项 (Note)

For research use only .