

产品名称: Cyclin B1 (Phospho Ser126) Rabbit Monoclonal Antibody
产品货号: AMRe21488



产品概述 (Summary)

产品名称 (Production Name)	Cyclin B1 (Phospho Ser126) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Phospho
修饰 (Modification)	Phosphorylated
同种型 (Isotype)	IgG,Kappa
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Protein A

免疫原信息 (Immunogen)

基因名 (Gene Name)	CCNB1
别名 (Alternative Names)	CCNB1;CCNB;G2/mitotic-specific cyclin-B1
基因 ID (Gene ID)	891.0
蛋白 ID (SwissProt ID)	P14635.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:500-1:2000,ICC/IF 1:200-1:1000,ELISA 1:5000-1:20000,IP 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:48kD;Observed MW:55kD


研究背景 (Background)

产品名称: Cyclin B1 (Phospho Ser126) Rabbit Monoclonal Antibody
产品货号: AMRe21488

Cell localization: Cytoplasm, Nuclear. The protein encoded by this gene is a regulatory protein involved in mitosis. The gene product complexes with p34(cdc2) to form the maturation-promoting factor (MPF). Two alternative transcripts have been found, a constitutively expressed transcript and a cell cycle-regulated transcript, that is expressed predominantly during G2/M phase. The different transcripts result from the use of alternate transcription initiation sites. [provided by RefSeq, Jul 2008],

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from HeLa

cells, using Cyclin B1 (p-Ser126) Rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG antibody was used to detect the antibody.

注意事项 (Note)

For research use only .