

产品名称: Iba1 Rabbit Monoclonal Antibody
产品货号: AMRe21499



产品概述 (Summary)

产品名称 (Production Name)	Iba1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG,Kappa
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Protein A

免疫原信息 (Immunogen)

基因名 (Gene Name)	AIF AIF 1;AIF-1;Aif1;AIF1 protein;AIF1_HUMAN;Allograft inflammatory factor 1;Allograft inflammatory factor 1 splice variant G;allograft inflammatory factor-1 splice variant Hara-1;balloon angioplasty responsive transcription;BART 1;G1;G1 putative splice variant of allograft inflammatory factor 1;IBA 1;IBA1;interferon gamma responsive transcript;Interferon responsive transcript 1;interferon responsive transcript factor 1;Ionized calcium binding adapter molecule 1;Ionized calcium-binding adapter molecule 1;ionized calcium-binding adapter molecule;IRT 1;IRT1;Microglia response factor;MRF1;Protein g1;
别名 (Alternative Names)	
基因 ID (Gene ID)	199.0
蛋白 ID (SwissProt ID)	P55008.

产品名称: Iba1 Rabbit Monoclonal Antibody
产品货号: AMRe21499



产品应用 (Application)


稀释比 (Dilution Ratio)	WB 1:1000-1:5000, IHC 1:2000-1:4000, ICC/IF 1:200-1:1000, ELISA 1:5000-1:20000, IP 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:17kD; Observed MW:17kD

研究背景 (Background)


Cell localization: Cytoplasm. This gene encodes a protein that binds actin and calcium. This gene is induced by cytokines and interferon and may promote macrophage activation and growth of vascular smooth muscle cells and T-lymphocytes. Polymorphisms in this gene may be associated with systemic sclerosis. Alternative splicing results in multiple transcript variants, but the full-length and coding nature of some of these variants is not certain. [provided by RefSeq, Jan 2016],

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of lysates from Mouse spleen cells, using Iba1 Rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG antibody was used to detect the antibody.




产品名称: Iba1 Rabbit Monoclonal Antibody
产品货号: AMRe21499



Immunohistochemical analysis of paraffin-embedded Human liver tissue. 1, primary Antibody was diluted at 1:200(4°C,overnight). 2, EDTA pH 9.0 was used for antibody retrieval(>98°C,20min). 3,Secondary antibody was diluted at 1:200(room temperature, 30min).



Immunohistochemical analysis of paraffin-embedded Mouse liver tissue. 1, primary Antibody was diluted at 1:200(4°C,overnight). 2, EDTA pH 9.0 was used for antibody retrieval(>98°C,20min). 3,Secondary antibody was diluted at 1:200(room temperature, 30min).



Immunohistochemical analysis of paraffin-embedded Rat liver tissue. 1, primary Antibody was diluted at 1:200(4°C,overnight). 2, EDTA pH 9.0 was used for antibody retrieval(>98°C,20min). 3,Secondary antibody was diluted at 1:200(room temperature, 30min).

注意事项 (Note)

For research use only .