

产品名称: Survivin Rabbit Monoclonal Antibody
产品货号: AMRe21501



产品概述 (Summary)

产品名称 (Production Name)	Survivin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ICC/IF,ELISA,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG,Kappa
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Protein A

免疫原信息 (Immunogen)

基因名 (Gene Name)	BIRC5
别名 (Alternative Names)	BIRC5;API4;IAP4;Baculoviral IAP repeat-containing protein 5;Apoptosis inhibitor 4;Apoptosis inhibitor survivin
基因 ID (Gene ID)	332.0
蛋白 ID (SwissProt ID)	O15392.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,ICC/IF 1:200-1:1000,ELISA 1:5000-1:20000,IP 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:16kD;Observed MW:16kD

研究背景 (Background)

Cell localization:Nucleus.This gene is a member of the inhibitor of apoptosis (IAP) gene family, which encode negative


产品名称: Survivin Rabbit Monoclonal Antibody
产品货号: AMRe21501



regulatory proteins that prevent apoptotic cell death. IAP family members usually contain multiple baculovirus IAP repeat (BIR) domains, but this gene encodes proteins with only a single BIR domain. The encoded proteins also lack a C-terminus RING finger domain. Gene expression is high during fetal development and in most tumors, yet low in adult tissues. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, Jun 2011],

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from Jurkat

cells, using Survivin Rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG antibody was used to detect the antibody.

注意事项 (Note)

For research use only .