

产品名称: Clathrin Heavy Chain Rabbit Monoclonal Antibody
产品货号: AMRe21547



产品概述 (Summary)

产品名称 (Production Name)	Clathrin Heavy Chain Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG,Kappa
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Protein A

免疫原信息 (Immunogen)

基因名 (Gene Name)	CLTC
别名 (Alternative Names)	CLH17 CLTCL2 KIAA0034
基因 ID (Gene ID)	1213.0
蛋白 ID (SwissProt ID)	Q00610.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:200,ICC/IF 1:200-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	Calculated MW:192kD;Observed MW:180kD

研究背景 (Background)

Cell localization:Cytoplasm.Clathrin is a major protein component of the cytoplasmic face of intracellular organelles, called coated vesicles and coated pits. These specialized organelles are involved in the intracellular trafficking of receptors and


产品名称: Clathrin Heavy Chain Rabbit Monoclonal Antibody
产品货号: AMRe21547



endocytosis of a variety of macromolecules. The basic subunit of the clathrin coat is composed of three heavy chains and three light chains. [provided by RefSeq, Jul 2008],

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from HeLa

cells, using Clathrin Heavy Chain 1 Rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG antibody was used to detect the antibody.

注意事项 (Note)

For research use only .