

产品名称: Nanog Rabbit Monoclonal Antibody
产品货号: AMRe21566



产品概述 (Summary)

产品名称 (Production Name)	Nanog Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG,Kappa
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein
纯化方式 (Purification)	Protein A

免疫原信息 (Immunogen)

基因名 (Gene Name)	NANOG
别名 (Alternative Names)	NANOG;Homeobox protein NANOG;Homeobox transcription factor Nanog;hNanog
基因 ID (Gene ID)	79923.0
蛋白 ID (SwissProt ID)	Q9H9S0.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:5000,IHC 1:200-1:1000,ICC/IF 1:200-1:1000,ELISA 1:5000-1:20000,IP 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW:42kD;Observed MW:42kD

研究背景 (Background)


产品名称: Nanog Rabbit Monoclonal Antibody
产品货号: AMRe21566



Cell localization:Nucleus.The protein encoded by this gene is a DNA binding homeobox transcription factor involved in embryonic stem (ES) cell proliferation, renewal, and pluripotency. The encoded protein can block ES cell differentiation and can also autorepress its own expression in differentiating cells. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2015],

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from NCCIT

cells, using Nanog Rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG antibody was used to detect the antibody.

注意事项 (Note)

For research use only .