

产品名称: RIP Rabbit Monoclonal Antibody  
产品货号: AMRe21588



## 产品概述 (Summary)

|                        |                                        |
|------------------------|----------------------------------------|
| 产品名称 (Production Name) | RIP Rabbit Monoclonal Antibody         |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody |
| 宿主 (Host)              | Rabbit                                 |
| 应用 (Application)       | WB,IHC,ICC/IF,ELISA,IP                 |
| 种属反应性 (Reactivity)     | Human                                  |

## 产品性能 (Performance)

|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                             |
| 修饰 (Modification)   | Unmodified                                                                               |
| 同种型 (Isotype)       | IgG,Kappa                                                                                |
| 克隆 (Clonality)      | Monoclonal                                                                               |
| 形式 (Form)           | Liquid                                                                                   |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |
| 储存溶液 (Buffer)       | PBS, 50% glycerol, 0.05% Proclin 300, 0.05% protective protein                           |
| 纯化方式 (Purification) | Protein A                                                                                |

## 免疫原信息 (Immunogen)

|                        |          |
|------------------------|----------|
| 基因名 (Gene Name)        | RIPK1    |
| 别名 (Alternative Names) | RIP RIP1 |
| 基因 ID (Gene ID)        | 8737.0   |
| 蛋白 ID (SwissProt ID)   | Q13546.  |

## 产品应用 (Application)

|                          |                                                                                          |
|--------------------------|------------------------------------------------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:1000-1:5000,IHC 1:200-1:1000,ICC/IF 1:200-1:1000,ELISA 1:5000-1:20000,IP 1:50-1:200 |
| 蛋白分子量 (Molecular Weight) | Calculated MW:76kD;Observed MW:76kD                                                      |

## 研究背景 (Background)

Cell localization:Cytoplasmic, Membranous.catalytic activity:ATP + a protein = ADP + a phosphoprotein.,function:Promotes


**产品名称: RIP Rabbit Monoclonal Antibody**  
**产品货号: AMRe21588**



apoptosis and activation of NF-kappa-B. Required for TNFRSF1A mediated activation of NF-kappa-B.,PTM:Autophosphorylated on serine and threonine residues.,PTM:Proteolytically cleaved by caspase-8 during TNF-induced apoptosis. Cleavage abolishes NF-kappa-B activation and enhances pro-apoptotic signaling through the TRADD-FADD interaction.,similarity:Belongs to the protein kinase superfamily. TKL Ser/Thr protein kinase family.,similarity:Contains 1 death domain.,similarity:Contains 1 protein kinase domain.,subunit:Binds to the death domain of TNFRSF6 and TRADD. Is recruited by TRADD to TNFRSF1A in a TNF-dependent process. Binds RIPK3, UBCE7IP1 isoform 3 (ZIN), EGFR, IKBKG, TRAF1, TRAF2 and TRAF3. Interacts with BNLF1. Interacts with SQSTM1 upon TNF-alpha stimulation. May interact with MAVS/IPS1.,


## 研究领域 (Research Area)

## 图片 (Image Data)




Western blot analysis of lysates from HEK293

cells, using RIP Rabbit mAb. The HRP-conjugated Goat anti-Rabbit IgG antibody was used to detect the antibody.



Immunohistochemical analysis of paraffin-embedded Rat kidney tissue. 1, primary Antibody was diluted at 1:200(4°C,overnight). 2, EDTA pH 9.0 was used for antibody retrieval(>98°C,20min). 3,Secondary antibody was diluted at 1:200(room temperature, 30min).

**产品名称: RIP Rabbit Monoclonal Antibody**  
**产品货号: AMRe21588**



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissue. 1, primary Antibody was diluted at 1:200(4°C,overnight). 2, EDTA pH 9.0 was used for antibody retrieval(>98°C,20min). 3,Secondary antibody was diluted at 1:200(room temperature, 30min).

#### **注意事项 (Note)**

For research use only .