

产品名称: Lamin B1 Rabbit Monoclonal Antibody
产品货号: AMRe83677



产品概述 (Summary)

产品名称 (Production Name)	Lamin B1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Lamin B1
别名 (Alternative Names)	LMN; ADLD; LMN2; LMNB; Lamin B1;;Lamin B1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P20700.A synthesized peptide derived from human Lamin B1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:5000-1:10000,IHC 1:100-1:500,ICC 1:100-1:500,IP 1:50-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 66 kDa ; Observed MW: 70 kDa

研究背景 (Background)

Lamins are intermediate filament proteins that assemble into a filamentous meshwork, and which constitute the major


产品名称: Lamin B1 Rabbit Monoclonal Antibody
产品货号: AMRe83677




components of the nuclear lamina, a fibrous layer on the nucleoplasmic side of the inner nuclear membrane.

研究领域 (Research Area)


图片 (Image Data)



All lanes use AB0054 Lamin B1 Antibody at 1:10000 dilution for 1 hour at room temperature.




All lanes use AB0054 Lamin B1 Antibody at 1:10000 dilution for 1 hour at room temperature.




All lanes use AB0054 Lamin B1 Antibody at 1:10000 dilution for 1 hour at room temperature.


产品名称: Lamin B1 Rabbit Monoclonal Antibody
产品货号: AMRe83677



Immunofluorescent analysis using the Antibody at 1:150 dilution.




Immunofluorescent analysis of HepG2 cells, using Lamin B1 Antibody.




Immunohistochemical analysis of paraffin-embedded Human non-Hodgkin's lymphoma, using the Antibody at 1:100 dilution.


产品名称: Lamin B1 Rabbit Monoclonal Antibody
产品货号: AMRe83677



Immunohistochemical analysis of paraffin-embedded Human glioblastoma, using the Antibody at 1:100 dilution.



Immunohistochemical analysis of paraffin-embedded Human prostate cancer, using the Antibody at 1:100 dilution.



Immunohistochemical analysis of paraffin-embedded rat liver, using Lamin B1 Antibody.

注意事项 (Note)

For research use only .