

产品名称: Lrp2/Megalin Rabbit Monoclonal Antibody
产品货号: AMRe83680



产品概述 (Summary)

产品名称 (Production Name)	Lrp2/Megalin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Lrp2/Megalin
别名 (Alternative Names)	Low-density lipoprotein receptor-related protein 2; LRP-2; Glycoprotein 330; gp330; Megalin; LRP2; DBS;;LRP2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P98164.A synthesized peptide derived from human LRP2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 522 kDa ; Observed MW: 280 kDa

研究背景 (Background)


产品名称: Lrp2/Megalin Rabbit Monoclonal Antibody
产品货号: AMRe83680




LRP2, also named as gp330 and Megalin, belongs to the LDLR family. It acts together with cubilin to mediate HDL endocytosis. LRP2 may participate in regulation of parathyroid-hormone and para-thyroid-hormone-related protein release. Defects in LRP2 are the cause of Donnai-Barrow syndrome (DBS). The antibody is specific to LRP2.

研究领域 (Research Area)

图片 (Image Data)



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Western blot analysis of Lrp2 expression in 293T cell lysate.

注意事项 (Note)

For research use only .