

产品名称: Catenin delta 1/CTNND1 Rabbit Monoclonal Antibody
产品货号: AMRe83683



产品概述 (Summary)

产品名称 (Production Name)	Catenin delta 1/CTNND1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Catenin delta 1/CTNND1
别名 (Alternative Names)	CAS; P120CAS; P120CTN; CTNND; p120;delta Catenin
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O60716.A synthesized peptide derived from human delta Catenin

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 108 kDa ; Observed MW: 105 kDa

产品名称: Catenin delta 1/CTNND1 Rabbit Monoclonal Antibody
产品货号: AMRe83683




研究背景 (Background)


Key regulator of cell-cell adhesion that associates with and regulates the cell adhesion properties of both C-, E- and N-cadherins, being critical for their surface stability. Beside cell-cell adhesion, regulates gene transcription through several transcription factors including ZBTB33/Kaiso2 and GLIS2, and the activity of Rho family GTPases and downstream cytoskeletal dynamics.

研究领域 (Research Area)

图片 (Image Data)




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.


产品名称: Catenin delta 1/CTNND1 Rabbit Monoclonal Antibody
产品货号: AMRe83683



Western blot analysis of delta1 Catenin expression in HeLa cell lysate.



Immunofluorescent analysis using the Antibody at 1:50 dilution.



Immunofluorescent analysis of MCF-7 cells, using delta 1 Catenin/p120 Catenin Antibody .

注意事项 (Note)

For research use only .