

产品名称: MEK1 Rabbit Monoclonal Antibody
产品货号: AMRe83685



产品概述 (Summary)

产品名称 (Production Name)	MEK1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat,Cow,Dog

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MEK1
别名 (Alternative Names)	MAPKK1; MEK1; MP2K1; PRKMK1; kinase MEK1; ERK activator kinase 1; MAP kinase kinase 1; MAP2K1; MAPK/ERK kinase 1; MAPKK 1;;MEK1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q02750.A synthesized peptide derived from human MEK1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 43 kDa ; Observed MW: 45 kDa

产品名称: MEK1 Rabbit Monoclonal Antibody
产品货号: AMRe83685




研究背景 (Background)


Microtubule plus-end-directed kinetochore motor which plays an important role in chromosome congression, microtubule-kinetochore conjugation and spindle assembly checkpoint activation. Drives chromosome congression (alignment of chromosomes at the spindle equator resulting in the formation of the metaphase plate) by mediating the lateral sliding of polar chromosomes along spindle microtubules towards the spindle equator and by aiding the establishment and maintenance of connections between kinetochores and spindle microtubules.

研究领域 (Research Area)

图片 (Image Data)




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.


产品名称: MEK1 Rabbit Monoclonal Antibody
产品货号: AMRe83685




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Immunoprecipitate (IP) analysis using the Antibody at 1:50 dilution. (wb at 1:6K dilution)




Immunofluorescent analysis using the Antibody at 1:50 dilution.




产品名称: MEK1 Rabbit Monoclonal Antibody
产品货号: AMRe83685



Immunofluorescent analysis using the Antibody at 1:50 dilution.



Immunofluorescent analysis using the Antibody at 1:150 dilution.



Immunofluorescent analysis using the Antibody at 1:150 dilution.

注意事项 (Note)

For research use only .