

产品名称: C10orf63 Rabbit Monoclonal Antibody
产品货号: AMRe83692



产品概述 (Summary)

产品名称 (Production Name)	C10orf63 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	C10orf63
别名 (Alternative Names)	C10orf63; ENKUR; Enkurin;; Enkurin
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q8TC29.A synthesized peptide derived from human Enkurin

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	29 kDa

研究背景 (Background)

Adapter that functions to localize a calcium-sensitive signal transduction machinery in sperm to a calcium-permeable ion


产品名称: C10orf63 Rabbit Monoclonal Antibody
产品货号: AMRe83692



channel. Microtubule inner protein (MIP) part of the dynein-decorated doublet microtubules (DMTs) in cilia axoneme, which is required for motile cilia beating.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of C10orf63 expression in Mouse testis lysate.

注意事项 (Note)

For research use only .