

产品名称: PTBP2 Rabbit Monoclonal Antibody
产品货号: AMRe83693



产品概述 (Summary)

产品名称 (Production Name)	PTBP2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PTBP2
别名 (Alternative Names)	brPTB; MIBP; nPTB; Polypyrimidine tract binding protein 2; PTB like protein; PTBP 2Ptbp2; Splicing regulator;;PTBP2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9UKA9.A synthesized peptide derived from human PTBP2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	57 kDa

研究背景 (Background)


产品名称: PTBP2 Rabbit Monoclonal Antibody
产品货号: AMRe83693




RNA-binding protein which binds to intronic polypyrimidine tracts and mediates negative regulation of exons splicing. May antagonize in a tissue-specific manner the ability of NOVA1 to activate exon selection. In addition to its function in pre-mRNA splicing, plays also a role in the regulation of translation.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.




Western blot analysis of PTBP2 in (1) Neuro-2a cell lysate; (2) HeLa cell lysate.


产品名称: PTBP2 Rabbit Monoclonal Antibody
产品货号: AMRe83693




Immunofluorescent analysis using the Antibody at 1:50 dilution.



Immunofluorescent analysis using the Antibody at 1:150 dilution.



Immunofluorescent analysis using the Antibody at 1:200 dilution.



产品名称: PTBP2 Rabbit Monoclonal Antibody
产品货号: AMRe83693



Immunohistochemical analysis of paraffin-embedded human colon, using PTBP2 Antibody.

注意事项 (Note)

For research use only .