

产品名称: MRRF Rabbit Monoclonal Antibody
产品货号: AMRe83697



产品概述 (Summary)

产品名称 (Production Name)	MRRF Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MRRF
别名 (Alternative Names)	MRFF; MTRRF; Ribosome recycling factor, mitochondrial; RRF;;MRRF
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q96E11.A synthesized peptide derived from human MRRF

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 29 kDa ; Observed MW: 24 kDa

研究背景 (Background)

Responsible for the disassembly of ribosomes from messenger RNA at the termination of mitochondrial protein


产品名称: MRRF Rabbit Monoclonal Antibody
产品货号: AMRe83697




biosynthesis. Acts in collaboration with GFM2. Promotes mitochondrial ribosome recycling by dissolution of intersubunit contacts.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Western blot analysis of MRRF expression in HeLa cell lysate.

产品名称: MRRF Rabbit Monoclonal Antibody
产品货号: AMRe83697



注意事项 (Note)

For research use only .