

产品名称: IL-33 Rabbit Monoclonal Antibody
产品货号: AMRe83700



产品概述 (Summary)

产品名称 (Production Name)	IL-33 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC
种属反应性 (Reactivity)	Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	IL-33
别名 (Alternative Names)	DVS27; IL 1F11 ;IL 33; IL1F11; IL33; Interleukin 1 family member 11; Interleukin 33; Interleukin 33 precursor; Interleukin33; NFEHEV; NFHEV; ;IL 33
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q8BVZ5(mouse)(human:O95760).A synthesized peptide derived from mouse IL 33

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 30 kDa ; Observed MW: 33 kDa

产品名称: IL-33 Rabbit Monoclonal Antibody
产品货号: AMRe83700




研究背景 (Background)

Cytokine that binds to and signals through the IL1RL1/ST2 receptor which in turn activates NF-kappa-B and MAPK signaling pathways in target cells. Involved in the maturation of Th2 cells inducing the secretion of T-helper type 2-associated cytokines.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of IL-33 expression in mouse lung cell lysate.

注意事项 (Note)

For research use only .