

产品名称: SNF5 Rabbit Monoclonal Antibody
产品货号: AMRe83702



产品概述 (Summary)

产品名称 (Production Name)	SNF5 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SNF5
别名 (Alternative Names)	SMARCB1; BAF47; hSNF5; INI1; RDT; RTPS1; Sfh1p; SMARCB1; SNF5 homolog; SNF5L1; Snr1; SWI/SNF comp;SMARCB1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q12824.A synthesized peptide derived from human SMARCB1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	44 kDa

研究背景 (Background)


产品名称: SNF5 Rabbit Monoclonal Antibody
产品货号: AMRe83702




The SWI-SNF complex is involved in the activation of transcription via the remodeling of nucleosome structure in an ATP-dependent manner. Brm (also designated SNF2 α) and Brg-1 (also designated SNF2 β) are the ATPase subunits of the mammalian SWI-SNF complex.

研究领域 (Research Area)

图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Western blot analysis of SNF5 in (1) HeLa cell lysate; (2) K562 cell lysate.

注意事项 (Note)

For research use only .