

产品名称: HIP1 Rabbit Monoclonal Antibody
产品货号: AMRe83704



产品概述 (Summary)

产品名称 (Production Name)	HIP1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, IP
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HIP1
别名 (Alternative Names)	HIP 1; hip1; HIP1/PDGFRB fusion gene; HIP1 ; Huntingtin interacting protein 1;;HIP 1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O00291.A synthesized peptide derived from human HIP 1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	116 kDa

研究背景 (Background)


产品名称: HIP1 Rabbit Monoclonal Antibody
产品货号: AMRe83704




Plays a role in clathrin-mediated endocytosis and trafficking .Involved in regulating AMPA receptor trafficking in the central nervous system in an NMDA-dependent manner.Regulates presynaptic nerve terminal activity.

研究领域 (Research Area)

图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Western blot analysis of HIP1 in Jurkat cell lysate.

注意事项 (Note)

For research use only .