

产品名称: GM130 Rabbit Monoclonal Antibody
产品货号: AMRe83706



产品概述 (Summary)

产品名称 (Production Name)	GM130 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat,Monkey,Cow,Dog

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GM130
别名 (Alternative Names)	GM130; Gm130 autoantigen; GOLGA 2; Golga2; Golgi autoantigen; Golgin 95;;Golgin subfamily A member 2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q08379.A synthesized peptide derived from human Golgin subfamily A member 2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 113 kDa ; Observed MW: 130 kDa

产品名称: GM130 Rabbit Monoclonal Antibody
产品货号: AMRe83706




研究背景 (Background)


Golgi auto-antigen; probably involved in maintaining cis-Golgi structure. The Golgi apparatus, which participates in glycosylation and transport of proteins and lipids in the secretory pathway, consists of a series of stacked cisternae (flattened membrane sacs). Interactions between the Golgi and microtubules are thought to be important for the reorganization of the Golgi after it fragments during mitosis. This gene encodes one of the golgins, a family of proteins localized to the Golgi.

研究领域 (Research Area)

图片 (Image Data)




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.


产品名称: GM130 Rabbit Monoclonal Antibody
产品货号: AMRe83706



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Immunofluorescent analysis using the Antibody at 1:200 dilution.



Immunohistochemical analysis of paraffin-embedded human cervix carcinoma, using GM130 Antibody.

注意事项 (Note)

For research use only .