

产品名称: Actin(Alpha-Actin) Rabbit Monoclonal Antibody
产品货号: AMRe83708



产品概述 (Summary)

产品名称 (Production Name)	Actin(Alpha-Actin) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Actin(Alpha-Actin)
别名 (Alternative Names)	ACTA, ACTA1, ACTS, Actin, alpha skeletal muscle, Alpha-actin 1;alpha Skeletal muscle actin
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P68133.A synthesized peptide derived from human alpha Skeletal muscle actin

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	42 kDa

产品名称: Actin(Alpha-Actin) Rabbit Monoclonal Antibody
产品货号: AMRe83708




研究背景 (Background)


Actins are highly conserved proteins that are involved in various types of cell motility and are ubiquitously expressed in all eukaryotic cells.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of Actin (Alpha-Actin) expression in (1) A431 cell lysate; (2) MCF-7 cell lysate.



Immunofluorescent analysis of HeLa cells, using Actin (Alpha-Actin) Antibody .



Immunohistochemical analysis of paraffin-embedded mouse muscle, using Actin (Alpha-Actin) Antibody.

产品名称: Actin(Alpha-Actin) Rabbit Monoclonal Antibody
产品货号: AMRe83708



注意事项 (Note)

For research use only .