

产品名称: TRAF6 Rabbit Monoclonal Antibody
产品货号: AMRe83710



产品概述 (Summary)

产品名称 (Production Name)	TRAF6 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TRAF6 TNF receptor-associated factor 6; E3 ubiquitin-protein ligase TRAF6;
别名 (Alternative Names)	Interleukin-1 signal transducer; RING finger protein 85; TRAF6; RNF85; TRAF 6; TRAF-6;;TRAF6
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9Y4K3.A synthesized peptide derived from human TRAF6

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 60 kDa ; Observed MW: 58 kDa

产品名称: TRAF6 Rabbit Monoclonal Antibody
产品货号: AMRe83710




研究背景 (Background)


E3 ubiquitin ligase that, together with UBE2N and UBE2V1, mediates the synthesis of 'Lys-63'-linked-polyubiquitin chains conjugated to proteins, such as ECSIT, IKKKG, IRAK1, AKT1 and AKT2. Also mediates ubiquitination of free/unanchored polyubiquitin chain that leads to MAP3K7 activation.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.




产品名称: TRAF6 Rabbit Monoclonal Antibody
产品货号: AMRe83710




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Immunohistochemical analysis of paraffin-embedded Rat ovary, using the Antibody at 1:100 dilution.




Immunohistochemical analysis of paraffin-embedded Rat stomach, using the Antibody at 1:100 dilution.




Immunohistochemical analysis of paraffin-embedded Human prostate cancer, using the Antibody at 1:100 dilution.


产品名称: TRAF6 Rabbit Monoclonal Antibody
产品货号: AMRe83710



Immunohistochemical analysis of paraffin-embedded Human astrocytoma, using the Antibody at 1:100 dilution.



Immunohistochemical analysis of paraffin-embedded Mouse skeletal muscle - gastrocnemius , using the Antibody at 1:100 dilution.



Immunohistochemical analysis of paraffin-embedded human colon carcinoma, using TRAF6 Antibody.

注意事项 (Note)

For research use only .