

产品名称: TRAF6 Rabbit Monoclonal Antibody
产品货号: AMRe83711



产品概述 (Summary)

产品名称 (Production Name)	TRAF6 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TRAF6 TNF receptor-associated factor 6; E3 ubiquitin-protein ligase TRAF6;
别名 (Alternative Names)	Interleukin-1 signal transducer; RING finger protein 85; TRAF6; RNF85; TRAF 6; TRAF-6;;TRAF6
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9Y4K3.A synthesized peptide derived from human TRAF6

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 60 kDa ; Observed MW: 58 kDa

产品名称: TRAF6 Rabbit Monoclonal Antibody
产品货号: AMRe83711




研究背景 (Background)


E3 ubiquitin ligase that, together with UBE2N and UBE2V1, mediates the synthesis of 'Lys-63'-linked-polyubiquitin chains conjugated to proteins, such as ECSIT, IKBKG, IRAK1, AKT1 and AKT2. Also mediates ubiquitination of free/unanchored polyubiquitin chain that leads to MAP3K7 activation.

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of TRAF6 expression in Jurkat cell lysate.



All lanes use the Antibody at 1:2k dilution for 1 hour at room temperature.



产品名称: TRAF6 Rabbit Monoclonal Antibody
产品货号: AMRe83711



Immunohistochemical analysis of paraffin-embedded human lung carcinoma, using TRAF6 Antibody.

注意事项 (Note)

For research use only .