

产品名称: Talin2 Rabbit Monoclonal Antibody
产品货号: AMRe83713



产品概述 (Summary)

产品名称 (Production Name)	Talin2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Talin2
别名 (Alternative Names)	DKFZp451B1011;ILWEQ;Talin-2;talin2;TLN 2;TLN2;Talin 2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9Y4G6.A synthesized peptide derived from human Talin 2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 272 kDa ; Observed MW: 270 kDa

研究背景 (Background)

As a major component of focal adhesion plaques that links integrin to the actin cytoskeleton, may play an important role in


产品名称: Talin2 Rabbit Monoclonal Antibody
产品货号: AMRe83713



cell adhesion. Recruits PIP5K1C to focal adhesion plaques and strongly activates its kinase activity.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of Talin 2 in expression HeLa cell lysate.

注意事项 (Note)

For research use only .