

产品名称: TMF Rabbit Monoclonal Antibody
产品货号: AMRe83717



产品概述 (Summary)

产品名称 (Production Name)	TMF Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TMF
别名 (Alternative Names)	ARA160; TMF; Tmf1;;TATA element modulatory factor
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P82094.A synthesized peptide derived from human TATA element modulatory factor

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 123 kDa ; Observed MW: 160 kDa

研究背景 (Background)


产品名称: TMF Rabbit Monoclonal Antibody
产品货号: AMRe83717



Potential coactivator of the androgen receptor. Mediates STAT3 degradation. May play critical roles in two RAB6-dependent retrograde transport processes: one from endosomes to the Golgi and the other from the Golgi to the ER. This protein binds the HIV-1 TATA element and inhibits transcriptional activation by the TATA-binding protein (TBP).

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of TMF expression in HeLa cell lysate.

注意事项 (Note)

For research use only .