

产品名称: TTF1 Rabbit Monoclonal Antibody
产品货号: AMRe83726



产品概述 (Summary)

产品名称 (Production Name)	TTF1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TTF1
别名 (Alternative Names)	BCH; BHC; NK-2; TEBP; TTF1; NKX2A; TITF1; TTF-1; NKX2.1;;Nkx 2.1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P43699.A synthesized peptide derived from human Nkx 2.1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	39 kDa

研究背景 (Background)


产品名称: TTF1 Rabbit Monoclonal Antibody
产品货号: AMRe83726




Transcription factor that binds and activates the promoter of thyroid specific genes such as thyroglobulin, thyroperoxidase, and thyrotropin receptor. Crucial in the maintenance of the thyroid differentiation phenotype. May play a role in lung development and surfactant homeostasis. Forms a regulatory loop with GRHL2 that coordinates lung epithelial cell morphogenesis and differentiation.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.




产品名称: TTF1 Rabbit Monoclonal Antibody
产品货号: AMRe83726



Western blot analysis of TTF1 expression in 293 cell lysate.



Immunofluorescent analysis of Hela cells, using TTF1 Antibody .



Immunohistochemical analysis of paraffin-embedded human thyroid, using TTF1 Antibody.

注意事项 (Note)

For research use only .