

产品名称: TCEB2 Rabbit Monoclonal Antibody
产品货号: AMRe83729



产品概述 (Summary)

产品名称 (Production Name)	TCEB2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TCEB2
别名 (Alternative Names)	EloB; Elongin B; SIII; SIII p18; Tceb2;;Elongin B
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q15370.A synthesized peptide derived from human Elongin B

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 13 kDa ; Observed MW: 16 kDa

研究背景 (Background)

SIII, also known as elongin, is a general transcription elongation factor that increases the RNA polymerase II transcription


产品名称: TCEB2 Rabbit Monoclonal Antibody
产品货号: AMRe83729




elongation past template-encoded arresting sites. Subunit A is transcriptionally active and its transcription activity is strongly enhanced by binding to the dimeric complex of the SIII regulatory subunits B and C (elongin BC complex)

研究领域 (Research Area)

图片 (Image Data)



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Western blot analysis of TCEB2 expression in Jurkat cell lysate.

注意事项 (Note)

For research use only .