

产品名称: RAIDD Rabbit Monoclonal Antibody
产品货号: AMRe83731



产品概述 (Summary)

产品名称 (Production Name)	RAIDD Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RAIDD
别名 (Alternative Names)	CRADD;MGC9163;RAIDD;Death adaptor molecule RAIDD;Death domain containing protein CRADD;;CRADD
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P78560.A synthesized peptide derived from human CRADD

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	23 kDa

产品名称: RAIDD Rabbit Monoclonal Antibody
产品货号: AMRe83731




研究背景 (Background)

Adapter protein that associates with PIDD1 and the caspase CASP2 to form the PIDDosome, a complex that activates CASP2 and triggers apoptosis.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of RAIDD expression in HeLa cell lysate.

注意事项 (Note)

For research use only .