

产品名称: CETP Rabbit Monoclonal Antibody
产品货号: AMRe83734



产品概述 (Summary)

产品名称 (Production Name)	CETP Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CETP
别名 (Alternative Names)	BPIFF; HDLCQ10; CETP; Lipid transfer protein I; Cholesteryl ester transfer protein;; CETP
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P11597.A synthesized peptide derived from human CETP

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	55 kDa

研究背景 (Background)


产品名称: CETP Rabbit Monoclonal Antibody
产品货号: AMRe83734



Involved in the transfer of neutral lipids, including cholesteryl ester and triglyceride, among lipoprotein particles. Allows the net movement of cholesteryl ester from high density lipoproteins/HDL to triglyceride-rich very low density lipoproteins/VLDL, and the equimolar transport of triglyceride from VLDL to HDL.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of CETP expression in human plasma lysate.

注意事项 (Note)

For research use only .