

产品名称: DNAJC15 Rabbit Monoclonal Antibody
产品货号: AMRe83735



产品概述 (Summary)

产品名称 (Production Name)	DNAJC15 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DNAJC15
别名 (Alternative Names)	HSD18; MCJ; Methylation controlled J protein; DNAJC15; DNAJD1; GIG22; DJC1;DNAJC15
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9Y5T4.A synthesized peptide derived from human DNAJC15

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	16 kDa

研究背景 (Background)


产品名称: DNAJC15 Rabbit Monoclonal Antibody
产品货号: AMRe83735




Negative regulator of the mitochondrial respiratory chain. Prevents mitochondrial hyperpolarization state and restricts mitochondrial generation of ATP. Acts as an import component of the TIM23 translocase complex. Stimulates the ATPase activity of HSPA9.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Western blot analysis of DNAJC15 expression in MCF-7 cell lysate.



Immunohistochemical analysis of paraffin-embedded human colon carcinoma, using DNAJC15 Antibody.

产品名称: DNAJC15 Rabbit Monoclonal Antibody
产品货号: AMRe83735



注意事项 (Note)

For research use only .