

产品名称: Bim Rabbit Monoclonal Antibody
产品货号: AMRe83736



产品概述 (Summary)

产品名称 (Production Name)	Bim Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Bim
别名 (Alternative Names)	BAM; BIM; BOD; BimL; BimS; BimEL; BIM-beta6; BIM-beta7; BIM-alpha6;BCL2L11;;BIM
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O43521.A synthesized peptide derived from human BIM

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	22 kDa

产品名称: Bim Rabbit Monoclonal Antibody
产品货号: AMRe83736




研究背景 (Background)


Induces apoptosis and anoikis. Isoform BimL is more potent than isoform BimEL. Isoform Bim-alpha1, isoform Bim-alpha2 and isoform Bim-alpha3 induce apoptosis, although less potent than isoform BimEL, isoform BimL and isoform BimS. Isoform Bim-gamma induces apoptosis. Isoform Bim-alpha3 induces apoptosis possibly through a caspase-mediated pathway. Isoform BimAC and isoform BimABC lack the ability to induce apoptosis.

研究领域 (Research Area)

图片 (Image Data)




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.




Western blot analysis of Bim expression in A431 cell lysate.


产品名称: Bim Rabbit Monoclonal Antibody
产品货号: AMRe83736




Immunofluorescent analysis using the Antibody at 1:150 dilution.



Immunofluorescent analysis of Raji cells, using Bim Antibody.




Immunohistochemical analysis of paraffin-embedded Human esophageal carcinoma, using the Antibody at 1:100 dilution.



产品名称: Bim Rabbit Monoclonal Antibody
产品货号: AMRe83736



Immunohistochemical analysis of paraffin-embedded human cervix cancer, using Bim Antibody.



Immunohistochemical analysis of paraffin-embedded Human prostate cancer, using the Antibody at 1:100 dilution.

注意事项 (Note)

For research use only .