

产品名称: RNF7 Rabbit Monoclonal Antibody
产品货号: AMRe83745



产品概述 (Summary)

产品名称 (Production Name)	RNF7 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RNF7
别名 (Alternative Names)	CKBBP1; Rbx2; RNF7; ROC2; SAG;;RNF7
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9UBF6.A synthesized peptide derived from human RNF7

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	13 kDa

研究背景 (Background)

Probable component of the SCF (SKP1-CUL1-F-box protein) E3 ubiquitin ligase complex which mediates the ubiquitination


产品名称: RNF7 Rabbit Monoclonal Antibody
产品货号: AMRe83745




and subsequent proteasomal degradation of target proteins involved in cell cycle progression, signal transduction and transcription. CRLs complexes and ARIH1 collaborate in tandem to mediate ubiquitination of target proteins, ARIH1 mediating addition of the first ubiquitin on CRLs targets.

研究领域 (Research Area)

图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Western blot analysis of RNF7 expression in HepG2 cell lysate.

注意事项 (Note)

For research use only .