

产品名称: CBR1 Rabbit Monoclonal Antibody
产品货号: AMRe83747



产品概述 (Summary)

产品名称 (Production Name)	CBR1 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CBR1 15 hydroxyprostaglandin dehydrogenase [NADP+]; Carbonyl reductase
别名 (Alternative Names)	[NADPH] 1; CBR1; CRN; NADPH dependent carbonyl reductase 1; Prostaglandin 9 ketoreductase; SDR21C1;;CBR1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P16152.A synthesized peptide derived from human CBR1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	30 kDa

产品名称: CBR1 Rabbit Monoclonal Antibody
产品货号: AMRe83747




研究背景 (Background)


NADPH-dependent reductase with broad substrate specificity. Catalyzes the reduction of a wide variety of carbonyl compounds including quinones, prostaglandins, menadione, plus various xenobiotics. Catalyzes the reduction of the antitumor anthracyclines doxorubicin and daunorubicin to the cardiotoxic compounds doxorubicinol and daunorubicinol.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



产品名称: CBR1 Rabbit Monoclonal Antibody
产品货号: AMRe83747



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

注意事项 (Note)

For research use only .