

**产品名称:** beta II Tubulin Rabbit Monoclonal Antibody  
**产品货号:** AMRe83748



## 产品概述 (Summary)

|                        |                                            |
|------------------------|--------------------------------------------|
| 产品名称 (Production Name) | beta II Tubulin Rabbit Monoclonal Antibody |
| 描述 (Description)       | Recombinant rabbit monoclonal antibody     |
| 宿主 (Host)              | Rabbit                                     |
| 应用 (Application)       | WB,IHC,ICC/IF,ICC,FC,IP                    |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat                            |

## 产品性能 (Performance)

|                     |                                                                                             |
|---------------------|---------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                |
| 修饰 (Modification)   | Unmodified                                                                                  |
| 同种型 (Isotype)       | IgG                                                                                         |
| 克隆 (Clonality)      | Monoclonal                                                                                  |
| 形式 (Form)           | Liquid                                                                                      |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.    |
| 储存溶液 (Buffer)       | Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol. |
| 纯化方式 (Purification) | Affinity Purification                                                                       |

## 免疫原信息 (Immunogen)

|                        |                                                                                                |
|------------------------|------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | beta II Tubulin                                                                                |
| 别名 (Alternative Names) | class II beta tubulin isotype; Beta2; Tubulin beta 2 chain; TUBB; TUBB2;;beta 4B/2A/2B Tubulin |
| 基因 ID (Gene ID)        |                                                                                                |
| 蛋白 ID (SwissProt ID)   | P68371/Q13885/Q9BVA1.A synthesized peptide derived from human beta 4B Tubulin                  |

## 产品应用 (Application)

|                          |                                                                                              |
|--------------------------|----------------------------------------------------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100,IP 1:20-1:50 |
| 蛋白分子量 (Molecular Weight) | 50 kDa                                                                                       |

产品名称: beta II Tubulin Rabbit Monoclonal Antibody  
产品货号: AMRe83748




## 研究背景 (Background)


Tubulin is the major constituent of microtubules. It binds two moles of GTP, one at an exchangeable site on the beta chain and one at a non-exchangeable site on the alpha chain.

## 研究领域 (Research Area)


## 图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.




**产品名称: beta II Tubulin Rabbit Monoclonal Antibody**  
**产品货号: AMRe83748**



---

All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Immunohistochemical analysis of paraffin-embedded mouse brain, using beta II Tubulin Antibody.

#### **注意事项 (Note)**

For research use only .