

产品名称: GAS2 Rabbit Monoclonal Antibody
产品货号: AMRe83750



产品概述 (Summary)

产品名称 (Production Name)	GAS2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GAS2
别名 (Alternative Names)	Gas2; Growth arrest specific 2; ;GAS2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O43903.A synthesized peptide derived from human GAS2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	35 kDa

研究背景 (Background)

May play a role in apoptosis by acting as a cell death substrate for caspases. Is cleaved during apoptosis and the cleaved


产品名称: GAS2 Rabbit Monoclonal Antibody
产品货号: AMRe83750




form induces dramatic rearrangements of the actin cytoskeleton and potent changes in the shape of the affected cells. May be involved in the membrane ruffling process (By similarity).

研究领域 (Research Area)

图片 (Image Data)



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Western blot analysis of GAS2 expression in Jurkat cell lysate.

注意事项 (Note)

For research use only .