

产品名称: CABP Rabbit Monoclonal Antibody
产品货号: AMRe83757



产品概述 (Summary)

产品名称 (Production Name)	CABP Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CABP
别名 (Alternative Names)	CaBP1; Calbrain; Calcium binding protein 1; Calcium binding protein 5; Caldendrin; HCALB BR;; Calcium binding protein 1
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9NZU7.A synthesized peptide derived from human Calcium binding protein 1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	Calculated MW: 40 kDa ; Observed MW: 45 kDa

产品名称: CABP Rabbit Monoclonal Antibody
产品货号: AMRe83757




研究背景 (Background)


Modulates calcium-dependent activity of inositol 1,4,5-triphosphate receptors. Inhibits agonist-induced intracellular calcium signaling. Enhances inactivation and does not support calcium-dependent facilitation of voltage-dependent P/Q-type calcium channels.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



产品名称: CABP Rabbit Monoclonal Antibody
产品货号: AMRe83757



Western blot analysis of CABP expression in U-87 MG cell lysate.



Immunofluorescent analysis using the Antibody at 1:50 dilution.

注意事项 (Note)

For research use only .