

产品名称: SHB Rabbit Monoclonal Antibody
产品货号: AMRe83758



产品概述 (Summary)

产品名称 (Production Name)	SHB Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SHB
别名 (Alternative Names)	Shb;;SHB
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q15464.A synthesized peptide derived from human SHB

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	55 kDa

研究背景 (Background)

Adapter protein which regulates several signal transduction cascades by linking activated receptors to downstream


产品名称: SHB Rabbit Monoclonal Antibody
产品货号: AMRe83758




signaling components. May play a role in angiogenesis by regulating FGFR1, VEGFR2 and PDGFR signaling. May also play a role in T-cell antigen receptor/TCR signaling, interleukin-2 signaling, apoptosis and neuronal cells differentiation by mediating basic-FGF and NGF-induced signaling cascades. May also regulate IRS1 and IRS2 signaling in insulin-producing cells.

研究领域 (Research Area)

图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Western blot analysis of SHB expression in K562 cell lysate.

注意事项 (Note)

For research use only .