

产品名称: Ceruloplasmin Rabbit Monoclonal Antibody
产品货号: AMRe83759



产品概述 (Summary)

产品名称 (Production Name)	Ceruloplasmin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Ceruloplasmin
别名 (Alternative Names)	Ceruloplasmin; CP2; Ferroxidase;;Ceruloplasmin
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P00450.A synthesized peptide derived from human Ceruloplasmin

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000
蛋白分子量 (Molecular Weight)	Calculated MW: 122 kDa ; Observed MW: 140 kDa

研究背景 (Background)

Multifunctional blue, copper-binding (6-7 atoms per molecule) glycoprotein. It has ferroxidase activity oxidizing Fe²⁺ to


产品名称: Ceruloplasmin Rabbit Monoclonal Antibody
产品货号: AMRe83759




Fe³⁺ without releasing radical oxygen species. It is involved in iron transport across the cell membrane.

研究领域 (Research Area)

图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Western blot analysis of Ceruloplasmin expression in Rat liver cell lysate.

注意事项 (Note)

For research use only .