

产品名称: Tbp7 Rabbit Monoclonal Antibody
产品货号: AMRe83760



产品概述 (Summary)

产品名称 (Production Name)	Tbp7 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Tbp7
别名 (Alternative Names)	MIP224; PSMC4; RPT3; Tat binding protein 7; TBP 7;; PSMC4
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P43686.A synthesized peptide derived from human PSMC4

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	47 kDa

研究背景 (Background)

The 26S protease is involved in the ATP-dependent degradation of ubiquitinated proteins. The regulatory (or ATPase)


产品名称: Tbp7 Rabbit Monoclonal Antibody
产品货号: AMRe83760




complex confers ATP dependency and substrate specificity to the 26S complex.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Western blot analysis of Tbp7 expression in HeLa cell lysate.

产品名称: Tbp7 Rabbit Monoclonal Antibody
产品货号: AMRe83760



注意事项 (Note)

For research use only .