

产品名称: Tripeptidyl peptidase 2 Rabbit Monoclonal Antibody
产品货号: AMRe83763



产品概述 (Summary)

产品名称 (Production Name)	Tripeptidyl peptidase 2 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Tripeptidyl peptidase 2
别名 (Alternative Names)	TPP II; TPP2; Tripeptidyl aminopeptidase; Tripeptidyl peptidase 2; Tripeptidyl peptidase II;;TPP2
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P29144.A synthesized peptide derived from human TPP2

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	138 kDa

产品名称: Tripeptidyl peptidase 2 Rabbit Monoclonal Antibody
产品货号: AMRe83763




研究背景 (Background)


Tripeptidyl peptidase II is a component of the proteolytic cascade acting downstream of the 26S proteasome in the ubiquitin-proteasome pathway. It may be able to complement the 26S proteasome function to some extent under conditions in which the latter is inhibited.

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of Tripeptidyl peptidase II expression in HepG2 cell lysate.



Immunofluorescent analysis using the Antibody at 1:50 dilution.

产品名称: Tripeptidyl peptidase 2 Rabbit Monoclonal Antibody
产品货号: AMRe83763



Immunofluorescent analysis using the Antibody at 1:50 dilution.

注意事项 (Note)

For research use only .