

产品名称: Acetyl-alpha Tubulin(Lys40) Rabbit Monoclonal Antibody
产品货号: AMRe83764

产品概述 (Summary)

产品名称 (Production Name)	Acetyl-alpha Tubulin(Lys40) Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC,IP
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Acetylated
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	alpha Tubulin(acetylK40)
别名 (Alternative Names)	Alpha-tubulin 1; TUBA1; TUBA1A ; Tubulin alpha 1 chain; Tubulin alpha;alpha 1B/4A Tubulin
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P68363/P68366.A synthesized peptide derived from human alpha 1B Tubulin

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	Calculated MW: 50 kDa ; Observed MW: 52 kDa

产品名称: Acetyl-alpha Tubulin(Lys40) Rabbit Monoclonal Antibody
产品货号: AMRe83764




研究背景 (Background)


Tubulin is the major constituent of microtubules. It binds two moles of GTP, one at an exchangeable site on the beta chain and one at a non-exchangeable site on the alpha chain.

研究领域 (Research Area)

图片 (Image Data)




All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.

产品名称: Acetyl-alpha Tubulin(Lys40) Rabbit Monoclonal
Antibody
产品货号: AMRe83764



Western blot analysis of alpha Tubulin (acetyl K40) expression in Mouse spleen cell lysate.

注意事项 (Note)

For research use only .