

产品名称: PUMA Rabbit Monoclonal Antibody
产品货号: AMRe83766



产品概述 (Summary)

产品名称 (Production Name)	PUMA Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ICC,FC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PUMA
别名 (Alternative Names)	BBC3; bcl2 binding component 3; jfy1; puma; p53 up-regulated modulator of apoptosis; PUMA alpha;;BBC3
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q9BXH1.A synthesized peptide derived from human BBC3

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,ICC/IF 1:50-1:200,ICC 1:50-1:200,FC 1:20-1:100
蛋白分子量 (Molecular Weight)	Calculated MW: 21 kDa ; Observed MW: 18 kDa

产品名称: PUMA Rabbit Monoclonal Antibody
产品货号: AMRe83766




研究背景 (Background)


Essential mediator of p53/TP53-dependent and p53/TP53-independent apoptosis. Promotes partial unfolding of BCL2L1 and dissociation of BCL2L1 from p53/TP53, releasing the bound p53/TP53 to induce apoptosis.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.




Western blot analysis of PUMA expression in K562 cell lysate.




Immunofluorescent analysis using the Antibody at 1:50 dilution.


产品名称: PUMA Rabbit Monoclonal Antibody
产品货号: AMRe83766




Immunofluorescent analysis using the Antibody at 1:150 dilution.



Immunofluorescent analysis using the Antibody at 1:200 dilution.



Immunofluorescent analysis of HeLa cells, using PUMA Antibody .



产品名称: PUMA Rabbit Monoclonal Antibody
产品货号: AMRe83766



Immunohistochemical analysis of paraffin-embedded human breast cancer, using PUMA Antibody.

注意事项 (Note)

For research use only .