

产品名称: HOXB9 Rabbit Monoclonal Antibody
产品货号: AMRe83767



产品概述 (Summary)

产品名称 (Production Name)	HOXB9 Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC, IP
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HOXB9
别名 (Alternative Names)	Homeobox protein Hox-B9; Hox-2E; Hox-2.5; HOXB9; HOX2E;;HOX B9
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P17482.A synthesized peptide derived from human HOX B9

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, ICC 1:50-1:200, IP 1:20-1:50
蛋白分子量 (Molecular Weight)	28 kDa

研究背景 (Background)

Recognizes the ribosomal RNA gene promoter and activates transcription mediated by RNA polymerase I (Pol I) through


产品名称: HOXB9 Rabbit Monoclonal Antibody
产品货号: AMRe83767



cooperative interactions with the transcription factor SL1/TIF-IB complex. It binds specifically to the upstream control element and can activate Pol I promoter escape.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of HOXB9 expression in MCF7 cell lysate.

注意事项 (Note)

For research use only .