

产品名称: Insulin Rabbit Monoclonal Antibody
产品货号: AMRe83768



产品概述 (Summary)

产品名称 (Production Name)	Insulin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ICC
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Insulin
别名 (Alternative Names)	ILPR; IRDN; IDDM2; MODY10; Insulin; Insulin; INS;; Insulin
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	P01308.A synthesized peptide derived from human Insulin

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:200, ICC/IF 1:50-1:200, ICC 1:50-1:200
蛋白分子量 (Molecular Weight)	12 kDa

研究背景 (Background)

Insulin decreases blood glucose concentration. It increases cell permeability to monosaccharides, amino acids and fatty


产品名称: Insulin Rabbit Monoclonal Antibody
产品货号: AMRe83768




acids. It accelerates glycolysis, the pentose phosphate cycle, and glycogen synthesis in liver.

研究领域 (Research Area)

图片 (Image Data)



Immunofluorescent analysis of BxPC-3 cells, using Insulin Antibody .



Immunohistochemical analysis of paraffin-embedded mouse pancreas, using Insulin Antibody.

注意事项 (Note)

For research use only .