

产品名称: MLKL Rabbit Monoclonal Antibody
产品货号: AMRe83772



产品概述 (Summary)

产品名称 (Production Name)	MLKL Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide, 0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	MLKL
别名 (Alternative Names)	Mixed lineage kinase domain-like protein; hMLKL; ;MLKL
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	Q8NB16.A synthesized peptide derived from human MLKL

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000, IHC 1:100-1:200
蛋白分子量 (Molecular Weight)	54 kDa

研究背景 (Background)

Pseudokinase that plays a key role in TNF-induced necroptosis, a programmed cell death process. Activated following


产品名称: MLKL Rabbit Monoclonal Antibody
产品货号: AMRe83772




phosphorylation by RIPK3, leading to homotrimerization, localization to the plasma membrane and execution of programmed necrosis characterized by calcium influx and plasma membrane damage.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.




All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Western blot analysis of MLKL expression in HUVEC cell lysate.

产品名称: MLKL Rabbit Monoclonal Antibody
产品货号: AMRe83772



Immunohistochemical analysis of paraffin-embedded human colon, using MLKL Antibody.

注意事项 (Note)

For research use only .