

产品名称: Survivin Rabbit Monoclonal Antibody
产品货号: AMRe83774



产品概述 (Summary)

产品名称 (Production Name)	Survivin Rabbit Monoclonal Antibody
描述 (Description)	Recombinant rabbit monoclonal antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,IP
种属反应性 (Reactivity)	Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Monoclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Purified antibody in PBS with 0.05% sodium azide,0.05% protective protein and 50% glycerol.
纯化方式 (Purification)	Affinity Purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	Survivin
别名 (Alternative Names)	API4;Apoptosis inhibitor survivin;BIRC 5;EPR 1;IAP4;SVV;TIAP;Survivin
基因 ID (Gene ID)	
蛋白 ID (SwissProt ID)	O70201(mouse)(human:O15392).A synthesized peptide derived from mouse Survivin

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:100-1:200,IP 1:20-1:50
蛋白分子量 (Molecular Weight)	16 kDa

研究背景 (Background)


产品名称: Survivin Rabbit Monoclonal Antibody
产品货号: AMRe83774




Multitasking protein that has dual roles in promoting cell proliferation and preventing apoptosis. Component of a chromosome passage protein complex (CPC) which is essential for chromosome alignment and segregation during mitosis and cytokinesis.

研究领域 (Research Area)


图片 (Image Data)



All lanes use the Antibody at 1:2K dilution for 1 hour at room temperature.



Western blot analysis of Survivin expression in Neuro-2a cell lysate.



Immunohistochemical analysis of paraffin-embedded mouse ovarian cancer, using Survivin Antibody.

产品名称: Survivin Rabbit Monoclonal Antibody
产品货号: AMRe83774



注意事项 (Note)

For research use only .